



# TECHNOPATH

CLINICAL DIAGNOSTICS

## Multichem<sup>®</sup> CSF

Providing Third-Party Cerebral Spinal Fluid QC  
in a Liquid Stable Format



### Features

- ◆ Human based proprietary formulation
- ◆ Mimics performance of patient samples.
- ◆ Targeted at key clinical decision points.
- ◆ Liquid for ease of use.

### Specifications

- ◆ 36 month closed vial stability at -20°C to -80°C.
- ◆ 30 day open vial stability at 2°C to 8°C.

Product	Description	Configuration	Part Code
Multichem CSF	Unassayed Bi-Level	2 × 6 × 2mL	CF100X
	Assayed Bi-Level	2 × 6 × 2mL	CF100A

# Multichem<sup>®</sup> CSF Analyte List

## Analytes

---

Glucose  
Lactate  
IgG  
Protein



Designed to complement and support TECHNOPATH's Multichem Quality Control (QC) product range.

---

IAMQC Software provides Laboratory Managers and Technologists with a range of QC software tools to analyse their QC results in real-time. IAMQC Software products allow users to automate, centralise and standardise QC processes in a laboratory setting. Our combination of software modules satisfy the varying levels of QC requirements in individual laboratories and are easily tailored to meet different QC management expectations.

For more information, visit [www.technopathcd.com](http://www.technopathcd.com)



**TECHNOPATH**  
CLINICAL DIAGNOSTICS

TECHNOPATH CLINICAL DIAGNOSTICS  
Technopath Life Sciences Park, Fort Henry, Ballina, Co. Tipperary, V94 FF1P, Ireland.  
For Multichem QC Support: [qcsupport@technopathcd.com](mailto:qcsupport@technopathcd.com) Website: [www.technopathcd.com](http://www.technopathcd.com)